



Twelve Sensors.  
Intelligently Positioned.

**CIRCA S-CATH™**  
Esophageal Temperature  
Monitoring System

**CIRCA S-CATH™ M\***  
Visible with 3-D Mapping

\*Pending 510(K) review. Not available in the United States.

FAST. ACCURATE. SENSITIVE.




**CIRCA**  
SCIENTIFIC

# RAPID, RESPONSIVE, CONTINUOUS MONITORING SOFTWARE

Continuous monitoring software is highly accurate in both hot and cold (down to 0°C) temperatures.<sup>1</sup>

- Four, user-selectable low and high temperature alarms
- Visual alarms for enhanced recognition
- Graphic and numeric temperature display
- Temperature log retains highest and lowest temperatures
- Conveniently record data for research



## STATIONARY PLACEMENT

Sensor placement ensures proximity to the point of treatment; no need to move the probe once placed.

- Radiopaque shaft provides a visual landmark of the esophagus
- Indicates esophageal width and orientation
- Facilitates reduced use of fluoroscopy



## CIRCA S-CATH™

Repositioned temperature sensors

Redesigned S-shape for expanded coverage

Pebax® coating electrically insulates - no bare metal

## CIRCA S-CATH™ M\*

Visible With 3-D Mapping

Minimal metal exposure on electrode

# CIRCA S-CATH™ TEMPERATURE PROBES

## Edge-to-Edge Coverage

During therapeutic procedures, esophageal temperatures can change quickly. The new and improved S-CATH provides faster, more accurate temperature detection.

- Soft, flexible self-expanding probe conforms to esophageal shape
- Proprietary sensor construction ensures rapid temperature transfer
- Delivers 240 data points per second; 12 temperature sensors update 20 times per second

# EDGE-TO-EDGE COVERAGE

The S-CATH, with its unique S-shaped design, deploys an array of 12 temperature sensors throughout the length and width of the esophagus, positioning sensors near the source of temperature changes. Independent research has shown that sensor distance has a great effect on temperatures recorded.<sup>2,3,4</sup>

■ Active Sensor


Basic Probe



3-Sensor Probe



S-CATH Probe





Average Esophageal Width = 18.9mm<sup>5</sup>

## BETTER VISUALIZATION

### S-CATH M\* is designed for 3-D Mapping Systems

- Four electrodes allow imaging on impedance-based 3-D mapping systems
- Electrodes centrally located to facilitate proper placement using 3-D mapping systems
- Twelve-sensor array provides temperature coverage without need to reposition



# IN VIVO DATA<sup>6</sup>

## Faster Detection

In an independent study of 198 applications in 10 patients, the S-CATH recognized an initial temperature rise of 0.2°C 17 seconds faster than a single sensor probe. (13.4±7.5 vs. 30.5±15.4 s; P, 0.001)

# 17 Seconds Faster



### INITIAL TEMPERATURE RISE:

**CIRCA** 17 Seconds Faster  
**S-CATH:** Giving you time to respond

## Multiple Sensors

In the same independent prospective study of 198 applications in 10 patients, a temperature rise of >2.0°C was recorded 40 times by the S-CATH. Single sensor probes missed 90% of those temperature rises.

# Single Sensor Missed 90%

### TEMPERATURE RISE >2.0°C RECORDED:

**CIRCA S-CATH:** 40  
**Single Sensor Probe:** 4

# BENCH DATA<sup>7</sup>

## Earlier Detection of a 2°C Rise


CIRCA S-CATH vs. single sensor 9F esophageal probe simultaneous submersion in warm water bath, representing optimal sensor positioning. Test conducted by CIRCA Scientific.

# 3X Faster

### TIME (SECONDS) TO DETECT A 2°C RISE

**CIRCA S-CATH:** 2.4  
**Single Sensor Probe:** 8.2

**Time to Reach 38°C**  
CIRCA S-CATH vs. Single Sensor Probe



Product Code	Description	
CS-1000	CIRCA Temperature Monitoring System™ (Touch Screen Display, Pole Mount Included)	CE 2797
CS-2006	CIRCA S-CATH™ Esophageal Temperature Probe (Single Use, 10Fr O.D., 10 units/Carton) International	CE 2797
CS-46EP*	CIRCA S-CATH™ M Esophageal Temperature Probe (Single Use, 10Fr. O.D., 10 units/Carton) International	CE 2797
CS-100*	CIRCA S-CATH™ M Mapping Interconnect Cable (Reuseable, 15 Foot Working Length)	CE
CS-101	CIRCA S-CATH™ Interconnect Cable (Reusable, 15 Foot Working Length)	CE
CS-1029	CIRCA Temperature Standard (Calibration)	CE
CS-1083	USB Data-Transfer Drive	CE

\*Pending 510(K) review. Not available in the United States.

**Indications for Use:** The CIRCA S-CATH and CIRCA S-CATH M Esophageal Temperature Probes are intended for continuous esophageal temperature monitoring during cardiac ablation procedures. The radiopaque probes are designed for placement in the esophagus. The CIRCA S-CATH M Esophageal Temperature Probe may be used to monitor electrophysiological signals. The CIRCA Temperature Monitor is indicated to display continuous temperature measurement (°C) from 12-sensor temperature probe for esophageal monitoring during cardiac ablation procedures.

1 Accuracy of the temperature sensors is  $\pm 0.3^{\circ}\text{C}$  within the rated output range of  $25^{\circ}\text{C}$  to  $45^{\circ}\text{C}$  and  $\pm 0.4^{\circ}\text{C}$  within the rated extended output range of  $0^{\circ}$  to  $24.9^{\circ}\text{C}$ .

2 Jose L. Merino, Martin Arceluz, Reina Delgado, Estela Falconi, Federico Cruz, Carlos C. Vasquez, Marta Ortega. Sensitivity and accuracy of Sensitherm/Esotherm oesophageal temperature probe: reply. *Europace* 2016;18:468-469.

3 Nakagawa H, Yamanashi WS, Pitha JV, Arruda M, Wang X, Ohtomo K, et al. Comparison of in vivo tissue temperature profile and lesion geometry for radiofrequency ablation with a saline-irrigated electrode versus temperature control in a canine thigh muscle preparation. *Circulation* 1995;91:2264-73.

4 Moreno J, Quintanilla JG, Molina-Morúa R, García-Torrent MJ, Angulo-Hernández MJ, Curiel Llamazares C, et al. Morphological and thermodynamic comparison of the lesions created by 4 open irrigated catheters in 2 experimental models. *J Cardiovasc Electrophysiol* 2014;25:1391-9. Medline.

5 Cury RC, Abbara S, Schmidt S, Malchano ZJ, Neuzil P, Weichet J, Ferencik M, et al. Relationship of the esophagus and aorta to the left atrium and pulmonary veins: Implications for catheter ablation of atrial fibrillation. *Heart Rhythm* 2005; 2:1317-1323.

6 Tschabrunn, CM, Silverstein J, Berzin T, Ellis E, Buxton AE, Josephson ME, Anter E. Comparison between single- and multi-sensor oesophageal temperature probes during atrial fibrillation ablation: thermodynamic characteristics, *Europace* 2015 doi:10.1093/europace/euu356.

7 Internal data. Test conducted by CIRCA Scientific.



#### Corporate Office

14 Inverness Drive East, Suite H-136

Englewood, CO 80112

www.CIRCASCIENTIFIC.com

Office: 1.303.951.8767 • Fax: 1.303.951.8769

info@circascientific.com

U.S. Patents 9,155,476 B2 and 9,668,655  
Other U.S. and foreign patents pending.

© 2019 CIRCA Scientific, Inc. All rights reserved.

CS-ART2072 Rev. 05