

## Product Description


The PeriCross™ Epicardial Access Kit is designed to provide safe and reliable access to the epicardium. The kit is supplied sterile and is designed to deliver a .018" guidewire to the pericardial space.

## System Components

- PeriCross™ Access Device
- PeriCross™ Tunneler Assembly
- Coaxial Introducer 5F x 18cm
- Guidewire – .018" OD x 80cm Long

## Ordering Information

The PeriCross™ Epicardial Access Kit  
Product Code: **PC-1013**



### Indications for Use

The PeriCross™ Epicardial Access Kit is intended to access the epicardial surface of the heart via a subxiphoid approach to facilitate electrophysiology studies in adult patients.

### RX Only Caution

Federal (U.S.A.) law restricts the use of this device for sale by or on the order of a physician.  
Patents: [www.circascientific.com](http://www.circascientific.com)

© 2026 CIRCA Scientific, Inc.

All rights reserved.

PeriCross is a trademark of CIRCA Scientific, Inc.† The third-party trademarks used herein are trademarks of their respective owners.