

S-CATH™ M Mapping Interconnect Cable

(CS-100)


PRODUCT SPECIFICATIONS

FEATURE	SPECIFICATIONS
Product Code	CS-100
Manufacturer	CIRCA Scientific, Inc.
Cable Length	4.6 m
Connectors	To: S-CATH M Probe Connector Type: Type C Linear Connector
	To: 3D Cardiac Mapping System Connector Type: 2mm Shrouded Pin Leads
Disposal	No special precautions required. Dispose of equipment per hospital policy.
Packaging	1 unit per polybag
Natural Rubber Latex Statement	Products and packaging are not made with natural rubber latex.
Sterilization	Not applicable
Warranty	2 years from date of sale
Electrical Safety Certification	CSA International

Connects S-CATH M Esophageal Temperature Probe to 3D Cardiac Mapping Systems

Other System Components:

- CS-21EP S-CATH M Esophageal Temperature Probe
- CS-1000 CIRCA Temperature Monitoring System
- CS-101 S-CATH Interconnect Cable
- CS-1029 CIRCA Temperature Standard (Calibration)



← To S-CATH M Esophageal Temperature Probe

To 3D Cardiac Mapping System →



CIRCA Scientific, Inc
Corporate Office
14 Inverness Drive East,
Suite H-136
Englewood, CO 80112
info@circascientific.com

Office: 1.303.951.8767
Fax: 1.303.951.8769